INSTRUCTION FOR USE

Resiment® Ready Mix with Fluoride bis-GMA Self-Curing Resin Cement

General Information

Resiment® Ready-Mix with Fluoride is a permanent bis-GMA auto-mixing, self-curing resin cement. Resiment® can be used for the final cementation of full-ceramic, porcelain fused to metal, composite, zirconia or full-cast alloy restorations, bridges, implant prosthesis, inlays and onlays, Maryland bridges, pins and posts. Resiment® is designed for easy use.

Advantages

- Resiment® Ready-Mix with Fluoride is an auto-mixing, auto-curing, filled, multi-purpose cement.
- Enjoy the convenience of an auto-mix delivery system.
- It can be used for final cementation of implant prostheses, bridges, crowns, inlays, Maryland bridges, cementation of pins, posts, pit and fissure sealant and periodontal splinting.
- Can be used for short clinical crowns.
- Can be used on dentin, enamel, composite, or metal.
- Is fluoride releasing.

Application

Initial Preparation of Dual Barrel Syringe

It is important to check for homogeneous and simultaneous expression of both base and catalyst components from their respective syringe orifice. Prior to initial use, remove cap from the unused syringe and slowly and evenly depress the plunger until equal portions of the base and catalyst components are expressed evenly and simultaneously from their respective openings, at the end of the double barrel. The system is now ready for use.

Secure the mixing tip onto the bayonet fixture of the Dual Barrel Syringe. Be sure to align the notch of the mixing tip with the groove on the attachment fixture of the syringe. Give the tip a quarter turn to the right to secure the tip in place. After use, dispose of the used mixing tip and be sure to replace with the original cap locked into place. If the cap is misplaced, an unused mixing tip can be left on the syringe until the next use. Do not leave the syringe barrel openings exposed for prolonged periods or during storage. Be careful not to cross-contaminate base and catalyst components at attachment openings.

Ceramic Restorations to Teeth

After try-in, clean by placing acid etched ceramic restoration (done previously per recommended manufacturer's indications and protocol*) in acetone or ethanol in ultrasonic cleaner for 1-2 minutes.

- Place silane coupling agent on intaglio surface of porcelain restoration.
- Acid etch (per manufacturer's indications and instructions) or clean smear layer of tooth with non-fluoridated prophy paste or flour of pumice.
- Prepare the surface of the tooth for bonding by using an enamel and dentin bonding system (per manufacturer's indications and instructions).
- Attach unused Resiment® Ready-Mix Tip to syringe and dispense mixture onto bonding surfaces to ensure all surfaces are covered.
- Fully seat the restoration and hold in place until Resiment® reaches initial set. (Approx. 1.5 minutes.)
- Remove excess cement IMMEDIATELY after initial set, using explorer, scaler, or brush. Do not wait until final set (3.5 minutes) to remove excess cement.

Zirconia Restorations to Teeth

Try in the restoration and verify margin fit, interproximal and occlusal contacts. Make final adjustments if necessary.

- Prepare the intaglio surface of the zirconia restoration by sandblasting with 50 micron aluminum oxide, clean surface with air and water.

- Prepare intaglio surface of restoration using an MDP containing priming agent.
- Thoroughly dry surfacing using air.
- Light cure primer inside of restoration. (Optional). Cover crown until ready to seat.
- Apply same primer to tooth surface. Dry with air and polymerize.
- Fill restoration wit Resiment® Ready-Mix with Fluoride.
- Fully seat restoration with pressure and verify margins for complete seating, have patient bite into cotton roll and hold in place until Resiment® reaches initial set. (Approx. 1.5 minutes.)
- Remove excess cement IMMEDIATELY after initial set, using explorer, scaler, or brush. Do not wait until final set (3.5 minutes) to remove excess cement.

Metal to Metal

Metal to metal restorations including implant prostheses, telescopic crown, or crown on post.

- Use a surface abrasive system (aluminum oxide sandblasting) on all external metal surfaces to be exposed to Resiment®. Intraoral metal should be roughened by a bur or a micro-etcher.
- If desired, apply a metal priming system of your choice to treat the metal exposed to Resiment®.
- Clean posts or abutments with non-fluoridated pumice prior to cementation to remove any debris from the surface.
- Attach unused Resiment® Ready-Mix with Fluoride tip to syringe and inject Resiment® into the restoration. Be sure not to trap any air bubbles.
- Fully seat restoration and hold in place until Resiment® reaches initial set. (Approx. 1.5 minutes)
- Remove excess cement IMMEDIATELY after initial set using explorer, scaler, or brush. Do not wait until final set (3.5 minutes) to remove excess cement.

Special Instructions

There always exists the remote possibility that a patient may have sensitivity towards Resiment® Ready-Mix with Fluoride. Should an allergic reaction occur, use of the material must be avoided and an alternative material, with different chemistry, should be employed.

Composition

Resiment® is a bis-GMA resin-based composite.

Storage

Do not store above 25°C/77°F

Keep away from children!

FOR PROFESSIONAL DENTAL USE ONLY.

Caution: Federal (USA) law prohibits this device to sale by or on the order of a dentist (or properly licensed practitioner).

Symbol	Used For	Symbol	Used For
	Batch code		Manufac-turer
	Do not reuse		Used by YYYY-MM-DD or YYYY-MM
	Sterilized using irradiation		Symbol for Non-Sterile

Symbol	Used For	Symbol	Used For
	Symbol for "Use by Prescription only"		Caution, consult accompanying documents
			Catalog number
			Do not use if package is damaged

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Sterngold Dental, LLC