# **TRU®/PUR®** Guided Surgery Kit

## Instructions for Use

23 Frank Mossberg Drive  
Attleboro, MA 02703  
+1.800.243.9942 or +1.508.226.5660

Order online at [www.sterngold.com](http://www.sterngold.com)

## US ENGLISH

### DESCRIPTION:

The Sterngold TRU® and PUR® Keyless Fully Guided Surgery Kit is a simple and predictable, surgical system that facilitates maximum ease of access by eliminating surgical keys. It is designed for dental implant placement of edentulous and partially edentulous jaws including single tooth replacement according to a case plan performed and/ or approved by the clinician. The kit is used exclusively to prepare the implant site and osteotomy for accurate placement of the full line of Sterngold TRU® & PUR® implants, both 3.5mm (Narrow Platform (NP)) and 4.5mm (Regular Platform (RP)) restorative platforms. Used in coordination with a previously designed and manufactured surgical guide, all drills and indicated surgical instruments are used through 5.04mm I.D. guide sleeves precisely placed into the guide based on the surgical planning. The guided drills are organized, and color coded according to implant diameter and are laser engraved by length to ensure ease of use during the drilling sequence. Other instruments to be used with the guide sleeves include both mounted and regular implant drivers, a guided flat drill for site preparation, NP and RP guided tissue punches, guided taps for each implant diameter, corresponding guided crestal bone drills for dense bone, and NP and RP guided latch drivers. Bone mills with NP and RP guide pins are included to modify the site prior to restoration. Additional instruments for use in the insertion and/or restorative processes include the thumb knob for use with the drivers, a 050 hex driver, a square to latch adaptor, a graduated torque wrench and guide fixation screws.

### INDICATIONS:

The Sterngold TRU® and PUR® Keyless Fully Guided Surgery Kit is to be used solely with Sterngold TRU® & PUR® implants. The kit is intended to be used for dental implant applications including oral rehabilitation of edentulous and partially dentate patients in the maxilla and mandible to support single unit, and multiple unit restorations. Implant retained restorations may consist of single crowns or bridges as well as complete or partial dentures. TRU® & PUR® implant is intended for delayed or immediate loading when adequate primary stability (35 Ncm) and appropriate occlusal loading is achieved clinically. Patient selection should be such that there are no compromising health conditions, disease, or medications that would interfere with the complete integration of the implants into the bone. Adequate bone quality and quantity should be assessed as well as available vertical opening to allow for kit instruments to be used.

### CONTRAINDICATIONS:

The following conditions would obviate the use of dental implants: inadequate bone quality and volume, infectable poor oral hygiene, acute or chronic infection, allergy to medical grade titanium, abuse of drugs or alcohol, cirrhosis, allergies to titanium, smoking, cardiopathy, neoplasia in action, systemic conditions that would impair healing, excessive occlusal parafunction, history of radiation or patient otherwise not suitable for long or complicated surgery, or the inability to consent to the procedure. The TRU® & PUR® Implant System is only intended for use with straight abutments.

TRU® & PUR® implant bodies are intended to be placed such that no angle correction is necessary.

### ⚠️ WARNINGS:

Ensure that each drill is used as directed by the surgical plan. Failure to do so may result in permanent damage and injury to the surgical site, nerves, or other vital anatomical structures. It is possible that an implant may fail to integrate with the surrounding bone, which would lead to loss of the implant and restoration supported by the implant. Potential causes for failure include inadequate surgical technique<sup>1</sup> not following the directions for use, poor case planning and/or patient selection which may include lack of bone quantity or quality, infection, poor patient oral hygiene or other co-morbid conditions. Temporary or permanent numbness (anesthesia), paresthesia, or dyesthesia are possible complications of implant surgery. Loosening or fracture of implant and restoration components may occur over time. Possible trauma may occur to implant site if incorrect selection of drill is used or used incorrectly, see below under Cautions.

### ⚠️ PATIENT PRECAUTIONS:

Following the surgical procedure and implant placement, the patient must carefully follow specific post-surgical directions including avoidance of prolonged strenuous physical activity following the procedure, consuming only soft foods, pain management and follow up as prescribed by the clinician.

by the surgical report and guide. According to surgical and sequence report, use appointed surgical components and drills to create osteotomy.

3. **Drilling Sequence:** The drilling sequence and specific size of each surgical component for each osteotomy will be indicated in the surgical and sequence report provided with the surgical guide. Note: TRU® and PUR® drills are not true-to-depth, the drill tip should not be included in determining the final depth of the osteotomy in coordination with indicator lines on each drill. Drilling speeds must be performed at speeds between 1500 - 2000 rpm for D1 and D2 bone and 1000 rpm for D3 and D4 w/copious irrigation refrigerated sterile saline (0.9%NaCl) This is a keyless surgical system; all components and drills will be advanced completely until they stop inside the surgical guide sleeves. The stops are built into the components and drills which are considered during the planning. All components and drills are color coded according to size and the surgical kit box. The general sequence of components and drills follow the following order:

- Tissue Punch for tooth and mucosa supported surgical guides not requiring flap technique allows for minimally invasive implant placement while preserving the adjacent tissue.
- Flat Drill, model and flatten the cortical bone at the implant site preparatory.
- 2.0mm x 6mm pilot drill. This will be the initial penetration of the cortical plate.
- Final osteotomy rill.
- Establish final width beginning with 2.0mm X 6mm then sequentially increase diameter until the final diameter is reached. Then proceed, increasing the length of the osteotomy using each sequential drill, until the final depth is reached. Continue advancing each drill until the stop on the drill contacts the guide sleeve.

Once the final diameter has been achieved, used the appointed drills within that diameter, continue with each consecutive depth drill in increasing order until the final appointed implant depth is reached.

7. Bone Tap, use the final appointed tap to create the threads for the implant.

8. Crestal Drill, use the appointed crestal drill for areas of high bone density to relieve the final diameter of the osteotomy in the crestal region.

- Implant Installation - TRU® and PUR® Implants** The TRU® & PUR® are self-tapping, double thread screw implants manufactured from Wrought Titanium (6AL-4V EL) as per ASTM F136. The implant surface is acid etched. The implants are manufactured with five body diameters 3.2 mm (only, PUR®), 3.5 mm, 4.3 mm, 5.0 mm, and 6.0 mm. The 3.2, 3.5, and 4.3 mm implants have a Narrow Platform (NP) prosthetic interface. The 5.0 and 6.0 mm implants have a Regular Platform (RP) prosthetic interface. The PUR® implants are substantially equivalent to the Zimvie Tapered Screw-Vent® 3.5 mm and the 4.5 mm platform implant system. The TRU® implants are substantially equivalent to the Nobel Connect® Conical Connection 3.5 mm and 4.5mm platform implant systems. They provide for non-ratlorial single and multiple tooth restorations in both the maxilla and mandible. They are packaged sterile, in a double blister which contains: the implant in a titanium tube and cover screw. Note: Refer to the Individual Instructions for use for specific warnings, cautions, and contraindications concerning TRU® and PUR® implants.
  - Using the guided driver of choice from the SGS® kit (Guided Latch, Guided Wrench, or Guided Mount Driver), install the TRU® & PUR® implants through surgical guide into the prepared osteotomy until the stop on the carrier/mount bottoms out. It is important to stop rotating the implant once the stop contacts the guide sleeve to avoid stripping the bone and any loss in primary stability. Using the graduated torque wrench, confirm primary stability of 35 N-cm, not to exceed 55 N-cm.

### ⚠️ CLEANING AND STERILIZATION:

Disinfection and sterilization procedures should conform to OSHA or local guidelines for blood borne pathogens.

### NON-STERILE RE-USABLE AND SINGLE USE COMPONENTS:

All non-sterile, re-usable tools and components must be cleaned and sterilized prior to use. These include: Titanium Forceps, Implant Driver, Implant Mount Drills, Bone Taps, Anchor Screws, Paralleling Pins, Pilot Drills, Lance Pilot Drills, Implant Driver Wrench, Implant Driver Latch, Drill Extender, 050° Latch Driver, 050° Hex Tools, Square to Latch Adaptor, Thumb Knobs, Graduated Torque Wrench, Torque Wrench, Ratchet Wrench, Lindemann Drill, Dense Bone Drills, Bone Mills, and Crestal Bone Drills. As well as all single use, non-sterile components including Abutments, Prosthetic Screws, Impression Copings, Implant Analogs, Scan Bodies, O-Rings, and Block out Spacers.

**CLEANING DISINFECTION PRIOR TO USE:** Devices are to be cleaned, disinfected, and sterilized by user prior to use. Use the following guidelines for cleaning/disinfecting for non-sterile re-usable components: Assemble multi-piece components, if applicable. Rinse with cool-to-lukewarm water for two-and-one-half minutes. For all parts place in an ultrasonic cleaner with an enzymatic detergent diluted with tap water per the manufacturer's guidelines. Sonicate for 10 minutes. Rinse with tap water for three minutes.

**INSPECTION AND FUNCTIONALITY:** To each use, special attention must be given to each instrument to confirm that it has been cleaned, visually inspected, and functionally tested to ensure it meets performance requirements as determined by the clinician. Once is has been determined that the instrument no longer meets these functional requirements, the instrument must be replaced.

All non-sterile, single-use components must be sterilized prior to use, with the exception of non-metal components, which should only be cleaned and disinfected.

**STERILIZATION:** Individual parts should be placed in appropriate autoclave for moist heat sterilization or in a dry heat pouch for dry heat sterilization. The following sterilization parameters (method, time, and temperature) are required to achieve a 10-6 sterility assurance level (SAL). Local or national specifications should be followed where steam sterilization requirements are stricter or more conservative than those listed in the table. Verify the calibration of your unit to ensure recommended temperatures are not being exceeded. To ensure autoclave is performing effectively, periodic use of biologic indicators should be considered.

Parts Individually produced	Cycle Type	Temperature	Exposure Time	Dry Time (only for kits)
	Gravity (steam)	121°C 250°F	40 minutes	n/a
Parts in a kit	Gravity (steam)	132°C 270°F	15 minutes	30 minutes
	Pre-vacuum (steam)	134°C 273°F	5 minutes	30 minutes

**⚠️ WARNING:** The use of non-sterile components may lead to infection of hard and/or soft tissue or may cause infectious disease or implant failure.

**⚠️ BREAKAGE:** Implant fractures can occur when applied loads exceed the normal functional design tolerances of the implant components. Potential overloading conditions may result from deficiencies in implant numbers, lengths and/or diameters to adequately support a restoration, excessive cantilever length, incomplete abutment seating, any abutment angulation or implant angle correction, occlusal interferences causing excessive lateral forces, patient parafunction (e.g., bruxing, clenching), improper denture manufacture procedures, inadequate denture fit, and physical trauma.

### CHANGES IN PERFORMANCE:

It is the responsibility of the clinician to instruct the patient on all appropriate contraindications, side effects, and precautions as well as the need to seek the services of a trained dental professional if there are any changes in the performance of the implant (e.g., looseness of the prosthesis, infection or exudate around the implant, pain, or any other unusual symptoms that the patient has not been told to expect).

### SHELF LIFE:

The product expiration date is indicated by the hourglass symbol on the product label, followed by the year and month of expiration.

**⚠️ CAUTION:** Federal (USA) law restricts this device to the sale by or on the order of a licensed clinician.

**HOW SUPPLIED:** Sterngold Dental implants are supplied sterile and should be handled with sterile titanium instruments. Implants are intended for single use only and should not be re-sterilized if contaminated.

### SOFTWARE:

These devices do not contain or utilize software.

### STORAGE:

Place devices in a dry place to prevent damage and/or deterioration.

### DISPOSAL:

Disposal of these components shall follow government regulations and environmental requirements.

**EMC AND ELECTRICAL SAFETY:** The TRU® & PUR® Implant System does not require EMC and Electrical Safety Evaluation.



### MR SAFETY:


The TRU® & PUR® Implant System has not been evaluated for safety and compatibility in the MR environment. The TRU® & PUR® Implant System has not been tested for heating or migration, or image artifact in the MR environment. The safety of TRU® & PUR® Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

**WARARRANTY STATEMENT:** Sterngold guarantees its products to be free from defects in material and/or workmanship. No other warranty is expressed or implied. This warranty applies to the original purchaser only. In the unlikely event of a defect, please follow the returned goods policy outlined in the product catalog.

### SYMBOLS:

The following table describes the symbols that may be printed on the packaging label. Please refer to the packaging label for the applicable symbols related to the product.

Symbol	Used For	Symbol	Used For
	Do Not re-use		Do not use if package is damaged
	Non-Sterile		Medical Device
	Sterilized using irradiation		UDI Unique Device Identifier
	Double sterile barrier system		Caution, consult accompanying documents
	Do not sterilize		Use by Date
	Use by Prescription Only		Part No.
	Consult Instructions for Use		Batch Code (lot #)

Symbol	Used For	Symbol	Used For
	Do Not re-use		Do not use if package is damaged
	Non-Sterile		Medical Device
	Sterilized using irradiation		UDI Unique Device Identifier
	Double sterile barrier system		Caution, consult accompanying documents
	Do not sterilize		Use by Date

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Sterngold Dental, LLC