SternVantage™ Monophase

ENGLISH

INSTRUCTION FOR USE

General

SternVantage™ Monophase is a medium viscosity thixotropic precision impression material in automix cartridges, based on addition curing vinyl silicones. SternVantage™ Monophase shows hydrophilic properties and thus improves the wettability and yields a better reproduction of subgingival areas.

SternVantage™ Monophase is thixotropic. Thus it stays in the tray but flows under pressure. Intra-oral tips are included for direct application on teeth.

Indication

One phase impression material for crowns, bridges, inlays, partial and complete dentures and denture repairs. Impression material for relining, functional impressions, impressions for prosthetic restorations and for impressions using the copper band technique.

Trays and Adhesive

Custom or non-perforated metal trays are recommended for best accuracy with SternVantage™ Monophase. Custom trays should not be used until 24 hours after manufacture.

In order to improve adhesion of SternVantage™ Monophase to the tray, SternVantage™ Adhesive can be applied as a thin even film. Let dry for 3 minutes.

Impression Taking

The material is loaded directly from the cartridge into the tray. Syringe material onto the prepared teeth until totally covered and insert the loaded tray no later than 1:30 minutes after start of mix into the mouth. Maintain position for at least 2:30 minutes.

Note: Times mentioned for 23°C (74°F). Higher temperatures shortens, lower temperatures prolong these times.

Casting and Galvanisation

The impression should be stored under normal room conditions (max. 25°C/ 77°F). The impression can be poured 30 minutes after removal from the mouth. There are no other time limitations.

All impressions can be copper- or silverplated.

Disinfection

Impressions can be disinfected with standard disinfection solutions.

Classification

SternVantage™ Monophase meets the requirements of ISO 4823 cat. A, type 2 (medium consistency).

Storage

Do not store above 25°C (77 °F). Do not use after expiration date.

Technical Data

Total working time (incl. mixing time)	1:30 min
Recommended time in mouth	2:30 min
Dimensional change after 24 h	≤ 0,5 %
Compression set	≤ 0,5 %
Strain in compression	≤ 3,5 %

Additional note: Do not use together with condensation curing silicone impression materials.

Packaging

SternVantage™ Monophase
2 cartridges @ 60 g, 12 mixing tips,
12 intra-oral tips

REF 220121

**Keep away from children!
For dental use only!**

Caution: Federal (USA) law prohibits this device to sale by or on the order of a dentist (or properly licensed practitioner).

Symbol	Used For	Symbol	Used For
	Do not reuse		Used by YYYY-MM-DD or YYYY-MM
	Batch code		Manufac-turer
	Sterilized using irradiation		Symbol for Non-Sterile
	Symbol for "Use by Prescription only"		Caution, consult accompanying documents
	Symbol for "European Conformity"		Catalog number
	Authorized representative in the European Community		

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Sterngold Dental, LLC