TempoLok® NE

ENGLISH

INSTRUCTION FOR USE

General Information

TempoLok® NE is an easy-to-use, eugenol-free based on zinc oxide, which does not affect the curing of resin-based filling and luting materials. TempoLok® NE flows easily, which allows perfect, accurate fitting of the temporary restoration. The optimized material properties of TempoLok® NE ensure reliable retention of the restoration while at the same time allowing easy removal. The excellent thermal insulation of TempoLok® NE prevents hot and cold sensitization. TempoLok® NE is mixed automatically in a double syringe and is applied directly to the restoration and into the cavity. Manual mixing of the two components is no longer required.

Indication

TempoLok® NE is used for temporary luting of crowns, inlays and onlays.

Handling of double syringe

The double syringe automatically dispenses and mixes the material, which saves time when applying the luting material directly into the restoration. Remove the transport cap before using a new double syringe and then attach the mixing tip. Following application, the mixing tip acts as a seal and can be left on the double syringe until the next use. The mixing tip should only be replaced by a new tip when the double syringe is used again.

Recommended use

Conditioning the preparation / restoration

Before cementation, clean the surface of the teeth to be restored (e.g. using pumice powder or water). Then rinse the tooth and dry with air. Avoid both excessive residual moisture and overdrying. The sensitivity of the prepared teeth can be reduced by using desensitizing agents (e.g. Fluoride Varnish). These agents should, however, be removed completely before permanent cementation. All residual cement also has to be carefully cleaned from the teeth when the temporary restoration is subsequently removed. Residual TempoLok® NE can be removed more easily, if the cement is cooled briefly using a cold spray.

Fitting the restoration

TempoLok® NE can be applied directly to the temporary restoration and into the cavity. Cement the temporary restoration onto the prepared tooth by applying light pressure. Then ask the patient to bite together slowly (use a cotton wool roll if necessary) and check the occlusion.

Removing excess material

TempoLok® NE is sufficiently cured after approx. 5 to 6 min. to allow removal of excess material. Begin removing excess cement proximally. Residual cement can be removed more easily, if the cement is cooled briefly using a cold spray. Then clean the interdental areas using dental floss.

Please note

Do not use TempoLok® NE in case of allergies to components of the material or the rare case of contact allergies are known. Avoid skin contact with the paste. On accidental contact, wash with soap and water. On contact with the eyes, rinse immediately with plenty of water and consult a doctor.

Working times

0 min	Application into the restoration or preparation.
0:30 min	Cementation of restoration onto prepared tooth.
1 min	End of the working time.
5 -6 min	Removal of excess material.
7 min	TempoLok® NE is fully cured.

Side effects

There are no known systemic side effects. In individual cases a hypersensitivity to components of the material cannot be ruled out.

Composition

Colophony, zinc oxide, natural resin, fatty acids, paraffin and in small amounts stabilizers and additives.

Classification

TempoLok® NE complies with the requirements of ISO 3107:2004 Dentistry – zinc oxide / eugenol and zinc oxide / noneugenol cements, Type 1, Class 1 (setting cement).

Storage

Store at normal room temperature (15-25 °C / 59-77 °F). Do not use after the expiration date.

Keep away from children! For dental use only!

Caution: Federal (USA) law prohibits this device to sale by or on the order of a dentist (or properly licensed practitioner).

Symbol	Used For	Symbol	Used For
	Do not reuse		Used by YYYY-MM-DD or YYYY-MM
	Batch code		Manufacturer
	Sterilized using irradiation		Symbol for Non-Sterile
	Symbol for "Use by Prescription only"		Caution, consult accompanying documents
	Catalog number		Do not use if package is damaged

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Sterngold Dental, LLC