prosthetic.line

Polyurock

DE	Gebrauchsanweisung	Deutsch	2
FR	Mode d'emploi	Français	3
EN	Instructions for Use	English	4
IT	Modo d'uso	Italiano	5
ES	Instrucciones de uso	Español	6

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Instructions for Use Polyurock

Composition

Polyurethane resin at 2 components, mixing proportions base : catalyst (5:1).

Physical and Mechanical properties at 25 °C (77 °F)

	Catalyst	Base	Mixture
Composition	Isocyanate	Polyol	
Color aspect	amber	white	white
Brookfield viscosity (mPa.s)	40	5	1.75
Density at 25 °C (g/cm ³)	1.20–1.24	1.62–1.68	1.55–1.61
Preparation time			~4 min 45
Final hardness (Shore D)			80–83
Modulus of elasticity (flexion) (Mpa)			3
Ultimate Tensile Strength (Mpa)			24
Compressive Yield Strength (Mpa)			62
Shrinkage (after 24 hours) (%)			0.07
Demoulding time (min.)			60–90

Instructions for Use

1. Mix colorant with base.
2. Thoroughly clean and dry the impression and box with wax or plasticine in the normal way.
3. Adapt Polyurock mixer on container of base and mix.
4. Fix the spout on the bottle of catalyst.
5. Pour out the desired amount of base to the blue line of the large goblet (9 or 10 for full models, 5 for half-models, small goblet for dies).
For small scales higher precision is obtained by weighing base to catalyst in a ratio 5:1.
6. Add catalyst to base up to corresponding red line of the goblet.
7. Mix for 30 seconds with a spatula.
8. Pour the mixture into the impression, use a vibrator to exclude air.
9. Possibility to put under pressure during 15 minutes at 2 bar.

Recommendations

Compatible with silicones, polyether (type Impregum), polyvinyle siloxane, alginate and hydrocolloid based impression material.
For polyethers (Impregum type), polyvinyl siloxanes, plastics and metals, use a commercial silicone spray for isolation.
Clean the impression with a freshly prepared solution containing 1 part household detergent to 3 parts distilled water and dry thoroughly.
No trace of moisture should remain.
Do not expose to steam more than 20 seconds (risk of deformation).
Avoid heating the die over 60 °C (140 °F).
Store in a dry place between 18 °C–20 °C (64 °F–68 °F).
Use rapidly after opening of container.
The model/die can be demoulded after 60–90 minutes.
Handle with care (wear gloves) and secure in container after use.
Clean with white spirit or acetone.

Caution

Do not eat, drink or smoke while using these products.
In case of contact with eyes rinse immediately with copious amounts of water and seek medical advice.
Both products, base and catalyst are irritating for the eyes, skin and respiratory system. They must be handled with care and secured in their container after use.
For further information see Material Safety Datasheet available from Cendres+Métaux SA, address here under.

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