**Bedienungsanleitung · Working instructions  
Mode d'emploi · Instrucciones de trabajo  
Istruzioni d'uso**

deutsch

english

français

español

italiano

# Dosper evo

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**Dosing and mixing unit for silicone masses**

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**Contents** (GB)

|       |                                                                                  |       |
|-------|----------------------------------------------------------------------------------|-------|
| 1.    | Technical Description                                                            | 22    |
| 2.    | Technical Data                                                                   | 22    |
| 3.    | Declaration of conformity                                                        | 23    |
| 4.    | Safety Advice                                                                    | 24    |
| 5.    | Initiation – Scope of Delivery                                                   | 25    |
| 6.    | Scheme of Functional Components                                                  | 26    |
| 7.    | Handling                                                                         | 27    |
| 7.1   | Connection of the Silicone Containers                                            | 27    |
| 7.1.1 | Connection of 2 x 1 kg (850 ml) Bottles                                          | 27    |
| 7.1.2 | Connection of 2 x 6 kg (5,1 L) and 2 x 25 kg (21,2 L)                            | 28    |
| 7.1.3 | Connection of the Material Reservoirs                                            | 28    |
| 7.2   | Connection of the Dosper evo to the Electric Circuit                             | 28    |
| 7.3   | Preliminaries                                                                    | 29    |
| 7.4   | Applying the Mixing Tips                                                         | 29    |
| 7.5   | Duplication                                                                      | 30    |
| 7.6   | Bottle-/Canister Change                                                          | 30    |
| 8.    | Deviations from the Regular Process                                              | 31    |
| 8.1   | Mixing Head only delivers one Component                                          | 31    |
| 8.2   | Proportion of the Components clearly deviates from a 1:1 Mixing                  | 32    |
| 8.3   | The Pump is in Operation, the LED's are Flashing, no Silicone is being delivered | 32    |
| 8.4   | The Machine shows no Function                                                    | 32    |
| 9.    | Optional Accessories                                                             | 33    |
| 10.   | Care and Maintenance                                                             | 33    |
| 11.   | Replacement of the Fuse                                                          | 34    |
| 12.   | General Advice                                                                   | 35    |
| 13.   | Guarantee                                                                        | 35    |
| 14.   | Spare Parts List                                                                 | 36    |
| 15.   | Exploded Drawing                                                                 | 37    |
| 16.   | Wiring Scheme                                                                    | 86-87 |

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
## 1. Technical Description

The Dosper evo is a mixing and dosing unit for all thin-flowing 1:1 duplicating silicones (addition-vulcanising or VPS) that are daily used in dental or orthodontic laboratories.

The machine delivers, doses, and mixes the silicones fully automatic and makes

the exactly required material quantity available at the push of a button.

Quick couplings allow the connection of all original packaging units (1kg, 6kg, 25 kg) as well as the optionally available clear-transparent material reservoirs

## 2. Technical Data

H x W x D (height x width x depth):

319 mm x 225 mm x 405 mm

Weight:

approx. 10 kg (without material)

Mains voltage:

220 – 230 V / 50-60 Hz 115 V / 60 Hz  
(Voltage / Frequency see type plate)

Power consumption:

approx. 280 W

Fuse:

1,6 A

Output:

150 ml / min.

Gear motor with overheat protection

Manual adjustability of delivered quantities

Illuminated work area

Optional:

Auto-stop before complete emptying  
(only when used with material reservoirs)

## 3. EC- Declaration of Conformity

### according to Low-Voltage Directive (2006/95/EG)

Herewith we declare that the machine described below in its conception and design and in the shape delivered by us is in accordance with the fundamental requirements for safety and health as prescribed in the applicable EC directives.

In the case of any change or modification of the machine, not authorized by us, this declaration becomes invalid.

### Denomination of the machine:

Dosper evo  
D-5050 / D-5050-A

### Type of machine:

Dosing and mixing unit for silicone masses

### Applicable harmonised standards:

2006/95/EG  
2004/108/EG

Low voltage guideline  
EMV-guideline

### Applicable harmonised standards 2006/95/EG

EN 61010-1:2001

### Low Voltage Directive

Safety requirements for electrical equipment for measurement, control and laboratory use - part 1: General requirements [IEC 61010-1:2001]

### 2004/108/EG

EN 61326-1:2006

### EMV-guideline

Electrical process measuring and control devices, laboratory units – EMV-requirements – part 1: General requirements (IEC 61326-1:2005)  
German Version 61326-1:2006  
Emission according to business sector and small enterprises

EN 55014-1S

Measuring of the interference voltage according to EN55014-1:2006 frequency range 150kHz – 30MHz  
Conducted emissions special evaluation of the interference

Dreve Dentamid GmbH  
Max-Planck-Straße 31 · 59423 Unna/Germany  
Tel.: +49 2303 8807-40 · Fax: +49 2303 8807-55  
E-Mail: info@dreve.de · www.dreve.com

## 4. Safety Advice

**Attention!** Please read these instructions carefully before connecting and putting into service. The reliability and the functionality of the machine can only be ensured, if the general governmental accident prevention regulations as well as the safety information in the manual are respected.

**1.** The machine may only be used according to the operating instructions on hand. We accept no liability for damages that are caused by improper handling or incorrect operation.

**2.** The Dosper must stand on an even undersurface which has to be sufficiently stable and resilient for the operating weight of 12-60 kg (depending on the used container sizes).

**3.** To avoid penetration of water (e.g. splashing water), the machine should be set up in a dry place.

**4.** The voltage mentioned on the type plate has to correspond with the voltage of the power source.

**5.** The machine shall only be connected to a socket with protective conductor. Never connect or disconnect the main plug with wet or moist hands.

**6.** Approved operators: the machine owner must make the user instructions accessible to the operators. He also has to make sure, that any approved operator has read and understood them. Only after that the operator is allowed to put the machine into operation.

**7.** The Dosper must be checked for its impeccable condition and work safety before use. If it is not in an impeccable condition, the machine shall not be used and has to be marked accordingly.

**8.** Do not put any objects into the machine.

**9.** Do not lock the pushbutton of the machine.

**10.** When not in use, the power supply has to be disconnected.

**11.** The signs and stickers have to be kept readable and shall not be removed.

**12.** Before cleaning and maintenance or replacement of spare parts disconnect the main plug.

**13.** The machine may only be opened and maintained by skilled personnel

**14.** Only use accessories and spare parts which have been authorized by the manufacturer. We cannot be held liable for any damage occurring from the use of unsuitable parts.

**15.** For safety reasons any unauthorized modification or change of the machine is prohibited.

**16.** The operating and maintenance instructions of this manual must abso-

lutely be followed. While working with the Dosper evo the general accident prevention regulations have to be kept to.

## 5. Initiation – Scope of Delivery

Before taking the machine into service, make sure that the mains voltage indicated on the type plate complies with the voltage of the power source.


Put the power cord into the appliance plug socket (3) and connect the Euro/UK/USA plug to the power source.

### Scope of delivery:

- 1 x Dosper evo
- 1 x power cord
- 1 x operating instructions in 5 languages
- 5 x mixing tips
- 1 x set bottle adapter

**Important:** Read the operating instructions carefully before initial startup. In case of uncertainties consult the responsible distributor or Dreve Dentamid GmbH directly.

## 6. Scheme of Functional Components

### No. Description

- |    |                                                                            |    |                   |
|----|----------------------------------------------------------------------------|----|-------------------|
| 1  | Main switch                                                                | 7  | Activation button |
| 2  | Fuse holder                                                                | 8  | LED lighting      |
| 3  | Mains connection                                                           | 9  | Setting screw     |
| 4a | Connector socket for autostop<br>(only if material reservoirs<br>are used) | 10 | Setting screw     |
| 4b | Connector socket for autostop<br>(only if material reservoirs<br>are used) |    |                   |
| 5a | Adapter for material container                                             |    |                   |
| 5b | Adapter for material container                                             |    |                   |
| 6  | Mixing head                                                                |    |                   |

## 7. Handling of the Dosper evo

Use high quality duplication silicones only. Unequal viscosities and distribution of filler materials in the individual components can lead to an abrasion in the pump housing and thus to a

deviation of the ideal 1:1 mixture (see point. 8.2)

We recommend the use of Dreve DUBLISIL-Silicones.

### 7.1 Connection of the Silicone Containers

With the help of adapters all commonly used material containers as well as the optionally available material reservoirs can be connected to the Dosper evo.

within the machine, causing major repairs which will have to be executed in our plant in Unna.

#### Important notice:

Please make sure that the component connections are never mixed up. On replacement of bottles or the refill of the reservoirs always ensure that the catalyst (standard colour white) as well as the base (coloured) are always connected to the same side of the machine.

#### Important notice:

The filling level of the containers must be controlled constantly! They shall never be emptied completely. Otherwise air might get into the system and the machine must be de-aerated in a complex procedure (see. point 8.1). A minimum level of about. 2-3 cm always has to remain in the containers.

Already through minimal blending of the components the silicone hardens

#### 7.1.1 Connection of 2 x 1 kg (850 ml) Bottles

Replace the original cap of the bottle by the adapter, which is included in the delivery, and fasten it hand-tight. Note that each adapter must absolutely be used for the same silicone component again when changing bottles in the future.

The bottles can now be connected to the machine with the help of the adapters. Turn the bottles upside down and

attach them to the respective counterparts of the adapters (5) on the machine. Press down until they noticeably snap into place.

To allow a problem-free flow of the silicone mass in the Dosper evo, and in order to avoid the development of a vacuum in the bottles, pierce tiny holes into their bottoms.

## 7.1.2 Connection of 2 x 6 kg (5,1 L) and 2 x 25 kg (21,2 L)

Replace the original caps of the canisters by the optionally available connection hoses (see point 9) and fasten them hand-tight.

Note that each connection hose must absolutely be used for the same silicone component again when changing canisters in the future.

The other side of the connection hose can now be connected to the machine with the help of the adapters. Turn the

canisters upside down and position them above the machine (a compatible support stand is also available (see point 9)). Press the connection hoses onto the adapters (5) until they noticeably snap into place.

To allow a problem-free flow of the silicone mass in the Dosper evo, and in order to avoid the development of a vacuum in the canisters, pierce tiny holes into their bottoms.

## 7.1.3 Connection of the Material Reservoirs

The optionally available material reservoirs (see point 9) are clear-transparent containers, marked with two different colours for distinction. They are connected to the machine with the help of the pre-assembled quick couplings. An electronically controlled float lever avoids the unwanted complete emptying of the reservoirs. To activate the system connect the plug-in contacts

mounted to the reservoirs with the respective sockets (4a/4b) of the machine.

Take off the lids of the reservoirs and fill in the respective silicone components. Each of the two reservoirs has a capacity of up to 1.5 L (corresponds to approx. 1.8 kg). Then close with the lids.

## 7.2 Connection of the Dosper evo to the Electric Circuit

Connect the provided power cord to the respective socket of the machine

(3), then connect to the electric circuit

## 7.3 Preliminaries

From the factory the Dosper evo is delivered with an inactive test liquid inside the system. This test liquid must completely be replaced by silicone before taking the machine into operation.

Turn the machine on at the main switch (1). The illuminated activation button (7) will light up.

Place an appropriate container underneath the mixing head and press the activation button (7).

The internal pump will press the test mass out of the machine and will replace it automatically by the connected silicone.

Keep the button (7) pressed until both silicone components are delivered evenly out of the mixing head (6).

**The Dosper evo is now ready for use**



## 7.4 Applying the Mixing Tips

The Dosper evo is equipped for the application of mixing tips with pink closure cap (D-4970-1).

There is a groove at the front of the mixing head (6). This groove has to

match with the triangular bulge of the mixing tip.

Snap the mixing tip into the guiding and apply to the mixing head (6) by a clockwise rotation of 90°.

## 7.5 Duplication

Place your duplicating flask under the mixing tip and press the activation button (7). The two LEDs (8) behind the mixing head will light up. On the one hand this is to control the operation of the pump, on the other hand the view in the work area is improved. As long as the silicone in the mixing tip is not cured, an arbitrary number of flasks can be filled one after the other. Leave the tip on the mixing head after the filling process, it will serve as a closure of the

mixing head. Next time the machine is going to be used just remove the old mixing tip by an anti-clockwise turn of 90° and replace it by a new one.

**Important:** Never activate the machine when a mixing tip with cured material is still locked to the mixing head!

In case the Dosper evo is not used for a longer period of time we recommend to turn the machine off at the main switch (1)

## 7.6 Bottle-/Canister Change

Release the quick couplings from the machine by pressing the locking button and pulling upwards. Screw the quick couplings to the new

bottle/canister as described above.

**Important:** Never swap the connection sides of the bottles/canisters (see chapter 8.1)


## 8. Deviations from the Regular Processes

### 8.1 Mixing Head only delivers one Component

**Cause:** As a bottle or a canister was completely emptied air has been sucked in by the pump. This might also happen when the minimum filling height falls below approx. 2 cm.

**Correction:** Replace the emptied silicone container and disassemble the

mixing head (see sequence of pictures) from the machine. Then put an appropriate container underneath the conveying hoses. Press the activation button until silicone comes out of both hoses equally. Reassemble the mixing head.

## 8.2 Proportion of the Components clearly deviates from a 1:1 Mixing.

**Cause:** The volume of one of the pump chambers has considerably increased due to abrasion – caused by quartz particles in the silicone.

**Correction:** Bring the filling height of the reservoirs (bottles/canisters) onto the same level and extract about 250 ml (corresponds to about 2-3 duplications). If the extraction volume after this try noticeably deviates in one of the components, the output can be adjusted to the required 1:1 dosage with the help of the setting screws (9/10)

**Important notice:** Exclusively change the output of the component, that shows a higher volume output (lower filling level)

Release the cover of the setting screw (9/10) that needs adjustment and carefully insert a slotted screwdriver into the opening until it noticeably snaps into the slot of the setting screw.

The volume output of this component can be reduced through anti-clockwise rotation.

Change the volume output only gradually. A ½ turn will already effect a volume reduction of about 4%. Repeat this process, if necessary, after you have checked the mixing ratio.

After adjustment re-lock the setting screw with its cover.

## 8.3 The Pump is in Operation, the LEDs are flashing, no Silicone is being delivered.

**Cause:** The silicone in the mixing tip is cured

**Correction:** Replace the cured mixing tip by a new one.

## 8.4 The Machine shows no Function

**Causes:**

- a) The machine is not connected to the power supply. Check the connection
- b) The main switch (1) is turned off. Turn the machine on.
- c) The sealing plugs respectively the connector plugs for the reservoirs (4a/4b) do not have contact. Pull the connectors from the machine, clean them, if necessary, and firmly insert them into their sockets.

d) The fuse is blown. Take the fuse out (see Chapter 11) and exchange, if necessary.

In case the Dosper evo still does not show any function, disconnect from the power supply and contact the responsible service technician or get in touch directly with the supplier Dreve.

## 9. Optional Accessories

**Material reservoirs:** D-5022

**Connection hose  
5,1 L canister:** D-49366

**Connection hose  
21,2 L canister:** D-49367

**Storage rack for  
5,1 L canister:** D-49364

**Storage rack for  
21,2 L canister:** D-49365

**Dublisil flask:** D-3805

**Mixing tips:**  
40 pieces D-4970-1  
100 pieces D-49701-1

## 10. Care and maintenance

Generally disconnect the unit from mains supply before executing care and maintenance.  
Ideally cleaning is to be effected dryly with a soft piece of cloth, if necessary


with a slightly humid sponge and a mild cleaning liquid. Water and cleaning liquids must not get into the unit.

## 11. Replacement of the Fuse

Disconnect the power cord from the current supply.

The fuse holder is integrated in the mains connection on the back of the machine.


For the control or the change of the fuse press the upper and the lower latch of the fuse holder simultaneously and take the holder out.

Pull the defective fuse out of the clamp holder by pulling gently. Replace with a new, identical fuse.

**Important notice:** Only use fuses with protection power (T 1, 6A) as stated on the back of the machine.

After the change insert the fuse holder again by pressing gently until the fixation locks noticeably.

## 12. General Advice

Repairs and maintenance shall only be executed by authorized and qualified personnel. In any case the electricity supply of the machine must be disconnected. The company DREVE is only responsible for the safety, reliability and performance of the device, if

- the machine is only used according to its purpose
- the machine is only operated by trained personnel

- the machine is used in compliance with this instructions manual
- enhancements, re-adjustments, changes or repairs are executed by persons authorized by the manufacturer Dreve

As we are constantly working on the advancement of our products, we reserve the right for technical modifications.

## 13. Guarantee

The right of the purchaser to bring actions arising out of defects shall in all cases be in lapse 2 years after the date of delivery. The guarantee obligation is only valid for material or production defects. It does not refer to natural wear and tear nor to damage occurring after the transfer of risk due to improper or unsuitable handling, excessive stressing, unsuitable operating material, and chemical, electrochemical or electrical influences of a nature not provided for according to the contract. On improper alterations or repairs or maintenance carried out by the purchaser or third parties, the guarantee obligations of the supplier and his liability for consequences caused thereby are rendered void.

### Liability

The right of the purchaser to bring actions arising out of defects shall in all cases be in lapse one year after the date of delivery. The warranty obligation does not refer to natural wear and tear or to damage occurring after the transfer of risk due to improper or unsuitable handling, excessive stressing, unsuitable operating material, and chemical, electrochemical or electrical influences of a nature not provided for according to the contract. On improper alterations, repairs or maintenance or breach of seals, carried out by the purchaser or third parties, our warranty obligations and liability for consequences caused thereby are rendered void.

## 14. Spare parts list


| No. | Denomination                                                          | Art.No.        |
|-----|-----------------------------------------------------------------------|----------------|
| 2   | Fuse                                                                  | 50048          |
| 3   | Combi element / mains connection                                      | 58514<br>58515 |
| 4   | Connector socket for auto stop (only if material reservoirs are used) | 58511          |
| 5a  | Adapter for material container                                        | 58508          |
| 5b  | Adapter for material container                                        | 58509          |
| 6   | Mixing head                                                           | D-6505045      |
| 7   | Activation button                                                     | 58513          |
| 8   | Set of LEDs                                                           | D-65058LED     |
| 10  | Setting screw / restrictor valve                                      | 58501          |
| 11  | Housing red                                                           | 58550rot       |
| 11a | Housing silver                                                        | 58550silber    |
| 11b | Housing special colour (indicate Ral-colour)                          | D-65058+RAL    |
| 12  | Bottom plate                                                          | 58555          |
| 13  | Appliance base foot (6 pcs.)                                          | D-6551079      |
| 14  | Motor & pump (indicate voltage)                                       | 50755          |
| 15  | Rubber buffer (4 pcs)                                                 | D-6558560      |
| 16  | Gasket ring for housing                                               | D-6558561      |
| 17  | Reducing nipple G3/8 a - G1/4i                                        | 58505          |
| 18  | Elbow connector 3/8 for hose 12/10                                    | 50264          |
| 19  | Polyamide hose 12x1                                                   | D-65028        |
| 20  | Straight screw joint G3/8                                             | 58502          |
| 21  | T-piece RP 3/8 I/I/I                                                  | 58500          |
| 22  | Double nipple 1/4 a - 3/8 a                                           | 50197          |
| 23  | Angular swivel joint F-WSI-1/4-MSv                                    | 58504          |
| 24  | Hollow bolt single 1/4                                                | 58506          |

|     |                                       |            |
|-----|---------------------------------------|------------|
| 25  | Straight screw joint-V.120614 6/4-1/4 | 50154      |
| 26  | Tecalan tube TR 6/4 mixing head       | D-6551058M |
| 27  | Tecalan tube TR 6/4 bypass            | D-6551068B |
| 28  | Angular screw joint G 3/8             | 58503      |
| 29  | Condensator 5µf 220-230V 50-60Hz EC   | 60808      |
| 29a | Condensator 20µf 115V 60Hz US         | 50796      |
| 29b | Condensator 25µf 100V 50/60Hz Japan   | 50797      |
| 30  | Electronics (indicate voltage)        | 58510      |
| 31  | Connecting terminal plate             | 51067      |
| 32  | Blind plug                            | 58512      |

### Without drawing No.

| Denomination                           | Art.No. |
|----------------------------------------|---------|
| Power cord Europe                      | 51021   |
| Power cord UK                          | 51213   |
|                                        | 51021   |
| Power cord USA/Japan/Korea             | D-65025 |
| PVC gasket 1/4                         | 50766   |
| FAG balls for mixing head              | 50271   |
| O-ring18x2 for bottle adapter          | 50769   |
| Pump chamber complete                  | D-65201 |
| Set of gear wheels including lock ring | D-65202 |
| Lid for pump flange complete           | D-65203 |
| Set of springs for mixing head         | D-65207 |
| Anti-interference capacitor            | 50799   |

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