



ERA® Implant Abutment

The smallest prosthetic head and the only with true vertical resiliency!

ERA® Abutments are manufactured for most popular screw and cylinder implants. Components and product numbers vary with the implant type. Refer to the ERA® ordering information on this page for compatibility. For ease in ordering, reference the group letter.

Call 800-243-9942 for assistance in ordering or go online to www.sterngold.com for an interactive ERA Implant Abutment Selection Chart. The direct link to the chart is: www.sterngold.com/Sterngold/Implants/ERA_Select.aspx



ERA® Implant Abutments fit the following implants:

SIZE	MANUFACTURER	GROUP
Sterngold-ImplaMed®		
3.3	Hex Cylinder	A
3.3	NP Self-tapping Hex Screw	W
3.3	NP Acid Etched	W
4.0	Hex Cylinder	A
3.75	Standard Hex Screw	A
3.75	Self-tapping Hex Screw	A
3.75	Self-tapping "SST" Hex Screw	A
4.0	Standard Hex Screw	A
4.0	Self-tapping Hex Screw	A
4.0	Self-tapping "SST" Hex Screw	A
5.0	RP "SST" Hex Screw	A
3.75	RP Acid Etched	A
4.0	RP Acid Etched	A
5.0	RP Acid Etched	A
5.0	WP Self-tapping Hex Screw	M
6.0	WP Self-tapping Hex Screw	M
3.3	Stern IC (4.8 head)	S
4.1	Stern IC (4.8 head)	S
Nobel Biocare Brånemark System®		
3.3	NP Nobel Speedy™ Groovy (ext.hex)	W
3.3	Fixture	A
3.3	NP Mk II	W
3.3	NP Mk III	W
3.75	Fixture	A
4.0	Fixture	A
5.0	Fixture (Old Version)	A
3.75	MkII, Self-tapping Fixture	A
4.0	MkII, Self-tapping Fixture	A
5.0	MkII, Self-tapping Fixture	M
5.5	MkII, Self-tapping Fixture	M
5.0	MkIV, Self-tapping Fixture	M
Nobel Biocare (Steri-Oss®)		
3.8	HL Cylinder	A
3.8	HL Threaded	A
4.5	HL Threaded	A
3.8	Cylindrical	D
3.8	Threaded	D
4.0	Steri-Oss	D
3.25	Cylindrical	L
3.5	Replace® Select & NobelReplace™	Z
4.3	Replace® Select & NobelReplace™	T
4.3	Nobel Active RP	B
5.0	Nobel Active RP	B
Interpore IMZ™		
3.3	Hex Cylinder	A
3.75	Self-tapping Threaded	A
4.0	Hex Cylinder	A
4.0	Self-tapping Threaded	A
4.25	Hex Cylinder	A
4.0	Cylinder	F
4.25	Cylinder	F
3.3	Cylinder	G
Straumann		
3.3	ITI TE™ (4.8 head)	S
3.3	ITI Std. & Std. Plus (4.8 head)	S
4.1	ITI TE™ (4.8 head)	S
4.1	ITI Std. & Std. Plus (4.8 head)	S
4.8	ITI Std. & Std. Plus (4.8 head)	S
Henry Schein®		
4.3	Camlog	AH
"O" Company		
3.25	Anti Rotational	E
4.0	Anti Rotational	D
Astra Tech		
3.5/4.0	Aqua	AK
4.5/5.0	Lilac	AJ

SIZE	MANUFACTURER	GROUP
BioHorizons®		
3.5	Internal	B
4.5	Internal	C
3.5	Single Stage	B
4.5	Single Stage	C
Osstem		
3.5	US II, II Plus	W
4.1	US II, III, II Plus, III Plus	A
	SI II, III (4.8 head)	S
3i Implant Innovations®		
3.25	External Hex Miniplant®	A
3.25	ICET™ Miniplant®	A
3.25	OSSEOTITE® Miniplant®	A
3.25	Internal Hex Miniplant®	E
3.3	Cylinder Miniplant®	A
3.3	External Hex Cylinder	A
3.75	ICET™ Self-tapping	A
3.75	OSSEOTITE®	A
3.75	Self-tapping Threaded	A
3.75	Standard Threaded	A
4.0	External Hex Cylinder	A
4.0	ICET™ Self-tapping	A
4.0	OSSEOTITE®	A
4.0	OSSEOTITE® Certain™	X
4.0	OSSEOTITE® NT Certain™	X
4.0	OSSEOTITE® CERTAIN PREVAIL	X
4.0	Standard Threaded	A
4.25	External Hex Cylinder	A
	TG OSSEOTITE® (4.8 Platform)	S
5.0	OSSEOTITE® Certain™	X
5.0	OSSEOTITE® NT Certain™	X
5.0	OSSEOTITE® CERTAIN PREVAIL	X
Zimmer (Paragon, Centerpulse)		
3.5	Bio-Vent® X™	A
3.75	Swede-Vent™ Conical Neck CST	A
3.75	Swede-Vent™ Standard	A
4.0	Swede-Vent™ Standard	A
4.0	Bio-Vent® X™	A
3.25	Micro-Vent® (3.5 head)	B
3.3	Screw-Vent® (3.5 head)	B
3.5	Bio-Vent® (3.5 head)	B
3.7	Screw-Vent® (3.5 head)	B
3.75	Screw-Vent® (3.5 head)	B
4.3	Core-Vent® (3.5 head)	B
4.25	Micro-Vent® (4.5 head)	C
4.5	Bio-Vent® (4.5 head)	C
4.7	Screw-Vent® (4.5 head)	C
5.3	Core-Vent® (4.5 head)	C
3.7	Tapered Swiss Plus™ (4.8 platform)	S
4.8	Tapered Swiss Plus™	S
4.1	Straight Swiss Plus™	S
4.8	Straight Swiss Plus™	S
Zimmer (Calcitek®, Centerpulse)		
3.25	Integral®	E
3.25	Omniloc®	E
3.75	ThreadLoc™	A
4.0	Integral®	D
4.0	Omniloc®	H
3.25	Spline®	J
3.75	Spline®	K
4.0	Spline®	K
Keystone (Lifecore)		
3.75	Restore® Self-tapping Screw	A
4.0	Restore® Self-tapping Screw	A
3.75	Restore® External Hex Screw	A
4.0	Restore® External Hex Screw	A
4.0	Restore® External Hex Cylinder	A
4.2	Sustain® External Hex Cylinder	A
3.75	Sustain® External Hex Screw	A

SIZE	MANUFACTURER	GROUP
4.0	Sustain® External Hex Screw	A
4.2	Sustain® External Hex MC Cylinder	A
4.0	Sustain® Internal Bevel	D
3.3	Stage-1™	S
4.0	Stage-1™	S
4.1	PrimaConnex®	AL
IMTEC Corporation®		
3.3	Universal Flare Cylinder	A
3.75	Universal Self-tapping	A
3.75	Universal Self-tapping Coated	A
4.0	Spike Cylinder	A
4.0	Universal Cylinder	A
Minimatic/Stryker		
3.3	External Hex Cylinder	A
3.75	External Hex Screw	A
4.0	External Hex Cylinder	A
4.0	External Hex Screw	A
4.75	External Hex Screw	A
5.0	External Hex Cylinder	A
OIC		
3.0	Osteo Standard ST	A
3.25	Osteo Standard ST	A
3.75	Osteo Standard ST	A
INNOVA		
4.1	ENDOPORE® Ext. Connection	A
4.0	ENTEGRAT™ Ext. Connection	A
Bud		
3.25	Bud Screwvent	A
3.75	Bud Screwvent	A
Biolok International		
4.5	Silhouette Screw	A
4.0	Micro-Lok Screw	A
4.0	Micro-Lok Cylinder	A
3.75	Micro-Lok Screw	A
3.3	Micro-Lok Cylinder	A
Implant Direct		
3.5	Legacy	B
4.5	Legacy	C
3.5	RePlant™	Z
4.3	RePlant™	T
3.7	ScrewPlant	B
4.7	ScrewPlant	C
Dentsply		
3.5	Ankylos	AE
4.5	Ankylos	AE
5.5	Ankylos	AE
3.8	FRIALIT® Plus	AC
3.8	XiVe® Plus	AC
3.8	XiVe® TG Plus	AC
4.5	FRIALIT® Plus	AD
4.5	XiVe® Plus	AD
4.5	XiVe® TG Plus	AD
MIS		
3.3	Internal Hex**	B
3.75	Internal Hex**	B
4.20	Internal Hex**	B
5.0	Internal Hex**	C
** these can be either the Biocom or Seven Implants		
ERA Micro or Standard Head		
ERA Micro Head only		
ERA Standard Head only		