

Now Serving...



SternVantage™ Quick Bite Vanilla

SternVantage™ Quick Bite Vanilla is a vanilla scented, automix thixotropic bite registration material based on addition curing polyvinyl siloxane.

Advantages

- Applied directly onto the teeth
- Allows a very fast bite registration (*1:45 minutes*)
- Hydrophilic – Compatible with moist surfaces
- Minimum resistance to the patient's jaw closure
- Produces consistently accurate centric relation records
- Enables the dentist and lab technician to seat models and dies easily into the bite registration
- Very rigid when set and can be milled and cut extremely well
- Can also be used as blockout material for attachments



SternVantage Quick Bite Vanilla 220131
*2 cartridges @ 60g each, 6 mixing tips,
6 contouration tips*

Mixing Tips 220162
50 pieces

Contouration Tips 220167
50 pieces

Dispensing Gun 220161
1:1 ratio



Order online at www.sterngold.com

Customer Service: 800-243-9942 • Toll Free Fax: 800-531-2685
23 Frank Mossberg Drive, P.O. Box 2967 • Attleboro, MA 02703-0967

**New
Product**