

Frequently Asked Questions

Aseptico AEU-1070SAU-70V

- **How long is the product warranty?**

On motors, two years against manufacturing defects. Handpieces are similarly covered for one year. Refer to the instruction manual for specific warranty policies.

- **Is the hand motor autoclavable?**

Yes, just make sure the red cap is attached to the black connector that goes into the console. Refer to the manual for maintenance instructions.

- **Is the tubing autoclavable?**

Yes, at least 90 times according to the manufacturer's instructions. Refer to the manual for maintenance instructions.

- **What handpieces does the unit support?**

There are settings for 20:1 (implant) and 1:1 (surgical) handpieces; however, any E-type handpiece may be used. Just set the unit to 1:1 and multiply or divide the display speed by the increaser or reduction ratio of the handpiece to know the bur speed.

- **Does the motor have handpiece calibration?**

Yes, for 20:1 implant handpieces.

- **Does the motor come in 220 V?**

The motor has automatic voltage sensing between 100 – 240 VAC, 50/60 Hz.

- **Where is the motor made and serviced?**

USA!

- **Whom do we call for service?**

Call Aseptico's Customer Service Team at 1-866-244-7095. Advanced Authorization is required before any unit is sent in for repair.