

Does your dental implant come with everything you need for a case?

TRU® and PUR® Implants

True Value. Pure Performance.



When it is time to restore a dental implant case, do you have all the necessary prosthetic components? With TRU® conical hex and PUR® internal hex connection implants, your restorative components are conveniently bundled with every implant, at a low flat-fee.

True value.

Convenience Bundles

Restorative components included with every implant*

Simple, Flat-Fee Pricing

PUR® or TRU® implant + restorative components

Performance

TRU® and PUR® implants are substantially equivalent⁽¹⁾ to the leading conical and screw-vent connection implants in the market. The precision manufacturing for TRU® and PUR® feature SGE surface treatment for expedited osseointegration.



⁽¹⁾The bone-level, internal connection TRU® implant is substantially equivalent to the Nobel Conical connection. Nobel is a registered trademark of Nobel Biocare™. The bone-level, internal connection PUR® implant is substantially equivalent to the Zimmer Tapered Screw-Vent®. Screw-Vent® is a registered trademark of Zimmer Biomet.

True Value

TRU® PUR®

Pure Performance

TRU®



PUR®



*Each implant is packaged with a cover screw

Pure performance.

TRU® conical and PUR® screw-vent implants are bone level, internal hex connection implants featuring a self-tapping, screw-type body design with a double lead thread pattern and micro threads at the top, ideal for achieving immediate stability.



TRU® Internal Conical Hex

- Simple, Predictable, Reliable
- Quality at True Value
- Provides Strong Initial Stability
- Diameter: 3.5mm, 4.3mm, 5.0 mm, 6.0mm⁽²⁾
- Lengths: 8mm, 10mm, 12mm, 14mm
- Mountless Delivery
- Comes in an All-in-1 Bundle
- Made in the USA

⁽²⁾ 5.0mm available in lengths of 8mm, 10mm, 12 mm and 6.0mm available in lengths of 8mm and 10mm



PUR® Internal Hex

- Simple, Predictable, Reliable
- Quality in Pure Performance
- Provides Strong Initial Stability
- Diameters: 3.2mm, 3.5mm, 4.3mm, 5.0mm, 6.0mm⁽²⁾
- Lengths: 8mm, 10mm, 12mm, 14mm
- Mountless Delivery
- Comes in an All-in-1 Bundle
- Made in the USA

Expect more.

Intuitively Designed Surgical Kit. An enhanced Surgical Kit featuring a comprehensive, easy drilling sequence, and color coded for handpiece and ratchet drivers. Crestal and dense bone drills featuring a matte finish for reduced glare with clear depth markings. Longer pins that can double up as depth gauges.


Add to Current Practice. TRU® and PUR® implants can be your implant of choice or can be added in to your current system as a value option.

Made in the USA. Manufactured under the strictest quality control processes at Sterngold Dental's US manufacturing facility in Massachusetts.

Surgical Kit



SGE Surface Treatment



Made in the USA



True Value

TRU® **PUR®**

Pure Performance

Convenience Bundles

All-In-1

Implant and restorative components in one affordable price. It's that simple. Sterngold offers analog, digital, or hybrid options.

Choose your implant*

TRU® or PUR®



+

Choose your restorative components

Select from our [analog](#) and/or [digital](#) prosthetics.

ANALOG



Healing Abutment



Open/Closed Tray Imp. Coping



Straight Abutment



Analog



Prosthetic Screw

OR

DIGITAL



Scan Body



Ti Base

OR

OR



Ti Base Prosthetic Screw

True Value

TRU® **PUR®**

Pure Performance

All for one bundled price

*Each implant is packaged with a cover screw

Convenience Bundles

The Practice

Lower your cost without sacrificing quality treatment for your patients.
Expect more from your dental implants.

TRU® or PUR® Implants + SURGICAL KIT

15 Implants + Kit

30 Implants + Kit

Each implant packaged with all restorative components of the All-in-1.



Optional Addition:

AEU-1070SAU-70V 70NCM Surgical Motor
& Mont Blanc 80 NCM 20:1 Contour Handpiece

