

TRU® and PUR® Implants

True Value. Pure Performance.



Convenience Bundles Order Form

All-in-1: One implant + Restorative Components

- (1) TRU®/PUR® Dental Implant*
- (1) TRU®/PUR® Healing Abutment
- (1) TRU®/PUR® Open/Closed Tray Impression Coping
- (1) TRU®/PUR® Analog
- (1) TRU®/PUR® Straight Final Abutment
- (1) TRU®/PUR® Prosthetic Screw

Digital option available: Replace straight abutment and impression coping with Ti Base & Scan Body



The Practice

BUNDLE IMPLANTS WITH SURGICAL KIT

OR GUIDED SURGERY KIT:

15 Implant Bundle + Kit

30 Implant Bundle + Kit



True Value

TRU® PUR®

Pure Performance

Select kit on page 7

*Each implant is packaged with a cover screw

Convenience Bundles

All-In-1

Implant and restorative components in one affordable price. It's that simple. Sterngold offers analog, digital, or hybrid options*.


Choose your implant*

TRU® or PUR®



Choose your restorative components

Select from our [analog](#) and/or [digital](#) prosthetics.



True Value

TRU® **PUR®**

Pure Performance

*Each implant is packaged with a cover screw



TRU® BUNDLE ORDER FORM

TRU® Narrow Platform (NP) Dental Implants (MAGENTA)
NP Prosthetic Connection = 3.5mm



*Each implant is packaged with a cover screw

3.5mm Body Diameter*	8mm	10mm	12mm	14mm	
Item #	901340	901341	901342	901343	TOTAL
Quantity					



TRU® NP Healing Abutments

Cuff Height	3mm	5mm	
Item #	901437	901439	TOTAL
Quantity			



TRU® NP Impression Copings

Item #	904015	TOTAL
Quantity		



TRU® NP Scan Body (Replaces Impression Copings)

Item #	904756	TOTAL
Quantity		



TRU NP Analog

Item #	904545	TOTAL
Quantity		



TRU® NP Straight Abutment

Cuff Height	1.5mm	3.0mm	
Item #	904003	904006	TOTAL
Quantity			



TRU® NP Ti Base (Replaces Straight Abutments)

Item #	904758	TOTAL
Quantity		



TRU® NP Prosthetic Screws

Item #	904656	TOTAL
Quantity		



Note: Total quantity of each TRU NP prosthetic should equal quantity of implants selected.
Ex: You select 30 implants, you must select 30 of each prosthetic in order to get bundle deal.



TRU® Regular Platform (RP) Dental Implants (YELLOW)

RP Prosthetic Connection = 4.5mm


*Each implant is packaged with a cover screw




4.3mm Body Diameter*	8mm	10mm	12mm	14mm	
Item #	901344	901345	901346	901347	TOTAL
Quantity					



5.0mm Body Diameter*					
Cuff Height	8mm	10mm	12mm		
Item #	901348	901349	901350	TOTAL	
Quantity					



6.0mm Body Diameter*	8mm	10mm		
Item #	901352	901353	TOTAL	
Quantity				



TRU® RP Healing Abutments

Cuff Height	3mm	5mm		
Item #	901438	901440	TOTAL	
Quantity				



TRU® RP Impression Copings

Item #	904017	TOTAL	
Quantity			



TRU® RP Scan Body (Replaces Impression Copings)

Item #	904757	TOTAL	
Quantity			



TRU® RP Analog

Item #	904544	TOTAL	
Quantity			



TRU® RP Straight Abutment

Cuff Height	1.5mm	3.0mm		
Item #	904005	904008	TOTAL	
Quantity				



TRU® RP Ti Base (Replaces Straight Abutments)

Item #	904759	TOTAL	
Quantity			



TRU® RP Prosthetic Screws

Item #	904657	TOTAL	
Quantity			



Note: Total quantity of each TRU RP prosthetic should equal quantity of implants selected.
Ex: You select 30 implants, you must select 30 of each prosthetic in order to get bundle deal.

PUR® BUNDLE ORDER FORM

PUR® Narrow Platform (NP) Dental Implants (SILVER)

NP Prosthetic Connection = 3.2mm on 3.2mm Implant, 3.5mm on 3.5mm/4.3mm implants

*Each implant is packaged with a cover screw

3.2mm Body Diameter*	8mm	10mm	12mm	14mm	
Item #	901382	901383	901390	901391	TOTAL
Quantity					



3.5mm Body Diameter*	8mm	10mm	12mm	14mm	
Cuff Height					
Item #	901356	901357	901358	901359	TOTAL
Quantity					



4.3mm Body Diameter*	8mm	10mm	12mm	14mm	
Item #	901360	901361	901362	901363	TOTAL
Quantity					



PUR® NP Healing Abutments

Cuff Height	3mm	5mm	
Item #	901446	901448	TOTAL
Quantity			



PUR® NP Impression Copings

Item #	904771	TOTAL
Quantity		



PUR® NP Scan Body (Replaces Impression Copings)

Item #	904790	TOTAL
Quantity		



PUR® NP Analog

Item #	904533	TOTAL
Quantity		



PUR® NP Straight Abutment

Cuff Height	1.5mm	3.0mm	
Item #	904773	904774	TOTAL
Quantity			



PUR® NP Ti Base (Replaces Straight Abutments)

Item #	904792	TOTAL
Quantity		



PUR® NP Prosthetic Screws

(Replaces NP Prosthetic Screw)

Item #	904653	Ti Base (904719)	TOTAL
Quantity			



Note: Total quantity of each PUR NP prosthetic should equal quantity of implants selected.
Ex: You select 30 implants, you must select 30 of each prosthetic in order to get bundle deal.

PUR®






PUR® Regular Platform (RP) Dental Implants (BLUE)

RP Prosthetic Connection = 4.5mm implants

*Each implant is packaged with a cover screw

5.0mm Body Diameter*	8mm	10mm	12mm	TOTAL	
Item #	901365	901366	901367		
Quantity					
6.0mm Body Diameter*	8mm	10mm	TOTAL		
Item #	901369	901370			
Quantity					

PUR® RP Healing Abutments

Cuff Height	3mm	5mm	TOTAL	
Item #	901447	901449		
Quantity				

PUR® RP Impression Copings

Item #	904772	TOTAL	
Quantity			

PUR® RP Scan Body (Replaces Impression Copings)

Item #	904791	TOTAL	
Quantity			

PUR® RP Analog

Item #	904534	TOTAL	
Quantity			

PUR® RP Straight Abutment

Cuff Height	1.5mm	3.0mm	TOTAL	
Item #	904775	904776		
Quantity				

PUR® RP Ti Base (Replaces Straight Abutments)

Item #	904793	TOTAL	
Quantity			

PUR® RP Prosthetic Screws

Item #	904653	(Replaces RP Prosthetic Screw) Ti Base (904719)	TOTAL	
Quantity				

Note: Total quantity of each PUR RP prosthetic should equal quantity of implants selected.
Ex: You select 30 implants, you must select 30 of each prosthetic in order to get bundle deal.

The Practice Bundle



TRU® or PUR® Implants + Surgical Kit or Guided Surgery Kit

15 Implants + Kit

30 Implants + Kit

Each implant packaged with all restorative components of the All-in-1.

OPTIONAL



AEU-1070SAU-70V 70NCM Surgical Motor & Mont
Blanc 80 NCM 20:1 Contour Handpiece

Complete Order Form

Name _____ Company Name _____

Address _____ Email _____

City _____ State _____ Zip _____

Office Telephone _____ Signature _____

Mobile Telephone _____

