


Stern IC (Internal Connection)

Features and Benefits

The **Stern IC** is an implant compatible with the **Straumann® ITI® Dental Implant**. It's single stage surgical procedure allows for reduced chair time and patient trauma.

- ◆ Mount body included
- ◆ 3.3 mm and 4.1 mm fixture diameters
- ◆ 1.8 mm machined collars
- ◆ 1.8 mm collar can also be used for 2 stage surgery
- ◆ Acid etched surface treatment
- ◆ Internal octagon
- ◆ Cover screw packaged separately to accommodate various restorative options
- ◆ Available lengths: 8 mm, 10 mm, 12 mm, and 14 mm



4.1 mm Reg. Neck, 1.8 mm Collar

8 mm	902721
10 mm	902722
12 mm	902723
14 mm	902724

3.3 mm Reg. Neck, 1.8 mm Collar

8 mm	902725
10 mm	902726
12 mm	902727
14 mm	902728

Healing Cap, Torx® Like

2 mm	901096
3 mm	901097

Analog

4 mm	904710
5.5 mm	904711
7 mm	904712

Cover Screw

Regular	901090
Torx® Like	901095

Waxing Sleeve

Kit	904497
Non-Engaging	904498
Engaging	904499

Solid Abutment

4 mm	904704
5.5 mm	904705
7 mm	904706

Impression Coping Reg. Neck

4 mm	904494
5.5 mm	904495
7 mm	904496



23 Frank Mossberg Drive • Attleboro, MA 02703
www.sterngold.com • 800-243-9942 • 508-226-5660

