

AEU-7000E



Aseptico's NEW Implant & Endodontic System

- Precision torque control
Endo: Auto-Stop-Reverse
Implant: Auto-Stop
- Upgradeable software
- Advanced calibration technology
- Programmable preset buttons



Model AEU-7000E-70V
Shown in Implant Mode
with Optional AHP-85MB
Handpiece



FEATURES	OPTIONAL ACCESSORIES
<ul style="list-style-type: none"> • Switch between Endo Mode and Implant Mode at the touch of a button • Powerful brushless 40,000 rpm autoclavable micromotor • Variable ratio selection: 1:5, 1:2, 1:1, 8:1 and 20:1 E-type handpieces • Advanced dynamometer calibration system accurately measures the properties of the micromotor and handpiece prior to use for ratios between: <ul style="list-style-type: none"> - Implant Mode 16:1 - 32:1 - Endo Mode 4:1 - 16:1 • Adjustable torque: <ul style="list-style-type: none"> - Implant Mode will reach torque up to 60 Ncm, which enables a single, high-efficiency 20:1 reduction handpiece (AHP-85MB) to be used for the complete implant procedure - Endo Mode will reach torque up to 1000 g-cm using 8:1 reduction handpiece • Selectable Auto-Stop function in Implant Mode and Auto-Stop-Reverse in Endo Mode when desired torque setting is reached • Contemporary design with easily cleaned surfaces 	 <ul style="list-style-type: none"> • AE-70V Variable Speed Foot Control (included)  <ul style="list-style-type: none"> • AE-7PM On/Off Foot Control (optional) • AHP-85MB 20:1 Reduction Contra Angle • AHP-85MB-C 20:1 Reduction Contra Angle • AHP-63 1:1 Standard Contra Angle • AHP-64 1:1 Standard Contra Angle • AHP-77C 1:2 Increaser Contra Angle • AHP-71TI 1:5 Increaser Contra Angle • AHP-88MN 8:1 Reduction Contra Angle • AHP-88MNP 8:1 Reduction Contra Angle
<p>FACTORY DEFAULT PRESETS</p> <p>Manual Mode: By default, the unit is always in manual mode. At any time, the user can adjust the torque, speed, irrigation flow, and other parameters, using the control panel keypad.</p> <p>Preset Mode: The system provides six memory presets for each Implant and Endo mode that can be used to quickly retrieve preferred settings for all operating parameters. Recalling these preset parameters saves time when preparing for different surgical/implant/endodontic procedures.</p> <p><i>NOTE: The default Preset Button settings will be restored when the unit's factory default settings are recalled or when the unit has been reprogrammed with new software.</i></p>	

REPLACEMENT PARTS			
Motor/Cable Assembly		Autoclavable Irrigation Tubing Set	Power Cord
AE-230-40	Brushless Micromotor	AE-23	840079
461561	Motor Holder	Peristaltic Pump Tubing Set (10pcs)	Memory Cards
461562	Motor Holder Attaching Bracket	AE-23-PUMP	MC-7000E
Footswitch		Irrigation Tubing Set for Bottles	Fuse
AE-7PM	On/Off Foot Control	AE-23-BOT	830040 SLO-BLO 1.60A 250V
Dynamometer Adapter		730576 <small>Replacement Disposable Bayonette for AE-23-BOT</small>	Manual for AEU-7000E-70V
461558		Irrigation Bag Hanger Rod	420763
		461541	



Address: 23 Frank Mossberg Drive • Attleboro, MA 02703
Toll-Free: 800-243-9942 • Phone: 508-226-5660
Web: www.sterngold.com • Fax: 508-226-5473