

ORA, Ball & O-Ring Compatibility Chart

IMPLANT DIAMETER	IMPLANT NAME	PROSTHETIC CONNECTION	GROUP	CUFF HEIGHT						
				0.4mm	1.0mm	1.25mm	2.0mm	3.0mm	4.0mm	5.0mm
BioHorizons®										
3.0	Internal	3.0 head	BH				904666	904667		
3.5	Internal	3.5 head	B		904340		904341	904342	904343	904344
4.0	Internal	4.5 head	C		904345		904346	904347	904348	904349
4.5	Internal	4.5 head	C		904345		904346	904347	904348	904349
3.5	Single Stage	3.5 head	B		904340		904341	904342	904343	904344
4.5	Single Stage	4.5 head	C		904345		904346	904347	904348	904349
4.0	External Hex (Maestro)	Regular	A			904335	904336	904337	904338	904339
Biolok International										
4.5	Silhouette Screw	Regular	A			904335	904336	904337	904338	904339
4.0	Micro-Lok Screw	Regular	A			904335	904336	904337	904338	904339
4.0	Micro-Lok Cylinder	Regular	A			904335	904336	904337	904338	904339
3.75	Micro-Lok Screw	Regular	A			904335	904336	904337	904338	904339
3.3	Micro-Lok Cylinder	Regular	A			904335	904336	904337	904338	904339
Biomet 3i										
3.25	External Hex Miniplant®	Regular	A			904335	904336	904337	904338	904339
3.25	ICE™ Miniplant®	Regular	A			904335	904336	904337	904338	904339
3.25	OSSEOTITE® Miniplant®	Regular	A			904335	904336	904337	904338	904339
3.3	Cylinder Miniplant®	Regular	A			904335	904336	904337	904338	904339
3.3	External Hex Cylinder	Regular	A			904335	904336	904337	904338	904339
3.75	ICE™ Self-tapping	Regular	A			904335	904336	904337	904338	904339
3.75	OSSEOTITE®	Regular	A			904335	904336	904337	904338	904339
3.75	Self-tapping Threaded	Regular	A			904335	904336	904337	904338	904339
3.75	Standard Threaded	Regular	A			904335	904336	904337	904338	904339
4.0	External Hex Cylinder	Regular	A			904335	904336	904337	904338	904339
4.0	ICE™ Self-tapping	Regular	A			904335	904336	904337	904338	904339
4.0	OSSEOTITE®	Regular	A			904335	904336	904337	904338	904339
4.0	Standard Threaded	Regular	A			904335	904336	904337	904338	904339
4.0	OSSEOTITE® Certain™		X		904390		904391	904392	904393	904394
4.0	OSSEOTITE® NT Certain™		X		904390		904391	904392	904393	904394
4.0	OSSEOTITE® Certain™ Prevail		X		904390		904391	904392	904393	904394
4.25	External Hex Cylinder	Regular	A			904335	904336	904337	904338	904339
3.25	TG OSSEOTITE®	4.8 head	S	904264	904265		904266	904267	904268	
4.0	TG OSSEOTITE®	4.8 head	S	904264	904265		904266	904267	904268	
5.0	TG OSSEOTITE®	4.8 head	S	904264	904265		904266	904267	904268	
5.0	OSSEOTITE® Certain™		X		904390		904391	904392	904393	904394
5.0	OSSEOTITE® NT Certain™		X		904390		904391	904392	904393	904394
5.0	OSSEOTITE® Certain Prevail		X		904390		904391	904392	904393	904394
Blue Sky Bio										
3.3	BIO One Stage	4.8 Head Regular	S	904264	904265		904266	904267	904268	
4.1	BIO One Stage	4.8 Head Regular	S	904264	904265		904266	904267	904268	
4.8	BIO One Stage	4.8 Head Regular	S	904264	904265		904266	904267	904268	
3.3	BIO One Stage	4.1 Head Narrow	S	904264	904265		904266	904267	904268	
4.1	BIO One Stage	4.1 Head Narrow	S	904264	904265		904266	904267	904268	
3.3	BIO Quattro	Narrow	BE		904550		904551	904552	904553	904554
4.1	BIO Quattro	Regular	BD		904469		904470	904471	904472	904474
4.8	BIO Quattro	Regular	BD		904469		904470	904471	904472	904474

IMPLANT DIAMETER	IMPLANT NAME	PROSTHETIC CONNECTION	GROUP	CUFF HEIGHT							
				0.4mm	1.0mm	1.25mm	2.0mm	3.0mm	4.0mm	5.0mm	
3.5	BIO Trilobe	Narrow	Z		904464		904465	904466	904467	904468	
4.3	BIO Trilobe	Regular	T		904374		904375	904376	904377	904378	
5.0	BIO Trilobe	Wide	AN		904670		904671	904672	904673	904674	
6.0	BIO Trilobe	Wide	AN		904670		904671	904672	904673	904674	
3.5	BIO Conus 12	3.5 / 4.0	AK		904690		904691	904692	904693	904694	
4.0	BIO Conus 12	3.5 / 4.0	AK		904690		904691	904692	904693	904694	
4.5	BIO Conus 12	4.5 / 5.0	AJ		904685		904686	904687	904688	904689	
5.0	BIO Conus 12	4.5 / 5.0	AJ		904685		904686	904687	904688	904689	
3.7	BIO Internal Hex	3.5	B		904340		904341	904342	904343	904344	
4.1	BIO Internal Hex	3.5	B		904340		904341	904342	904343	904344	
4.7	BIO Internal Hex	3.5	B		904340		904341	904342	904343	904344	
4.7	BIO Internal Hex	4.5	C		904345		904346	904347	904348	904349	
5.7	BIO Internal Hex	4.5	C		904345		904346	904347	904348	904349	
7.0	BIO Molar	4.5	C		904345		904346	904347	904348	904349	
3.5	BIO Max	Narrow	AY		904680		904681	904682	904683	904684	
4.3	BIO Max	Narrow	AY		904680		904681	904682	904683	904684	
5.0	BIO Max	Narrow	AY		904680		904681	904682	904683	904684	
Bud											
3.25	Bud Screwvent	Regular	A				904335	904336	904337	904338	904339
3.75	Bud Screwvent	Regular	A				904335	904336	904337	904338	904339
Dentsply											
3.5	Ankylos		AE		904555		904556	904557	904558	904559	
4.5	Ankylos		AE		904555		904556	904557	904558	904559	
5.5	Ankylos		AE		904555		904556	904557	904558	904559	
3.5	Astra	Aqua	AK		904690		904691	904692	904693	904694	
4.0	Astra	Aqua	AK		904690		904691	904692	904693	904694	
4.5	Astra	Lilac	AJ		904685		904686	904687	904688	904689	
5.0	Astra	Lilac	AJ		904685		904686	904687	904688	904689	
3.5 S	Astra	OsseoSpeed™ TX	AK		904690		904691	904692	904693	904694	
4.0 S	Astra	OsseoSpeed™ TX	AK		904690		904691	904692	904693	904694	
4.5	Astra	OsseoSpeed™ TX	AJ		904685		904686	904687	904688	904689	
5.0	Astra	OsseoSpeed™ TX	AJ		904685		904686	904687	904688	904689	
5.0 S	Astra	OsseoSpeed™ TX	AJ		904685		904686	904687	904688	904689	
Hiossen (Osstem)											
4.1	US II, III, II Plus, III Plus	Regular	A				904335	904336	904337	904338	904339
4.0	SI II, III	4.8 head	S	904264	904265		904266	904267	904268		
4.5	SI II, III	4.8 head	S	904264	904265		904266	904267	904268		
5.0	SI II, III	4.8 head	S	904264	904265		904266	904267	904268		
Implant Direct											
3.5	Legacy	3.5 head	B		904340		904341	904342	904343	904344	
4.5	Legacy	4.5 head	C		904345		904346	904347	904348	904349	
5.7	Legacy	5.7 head	BF		904695		904696	904697	904698	904699	
3.5	RePlant™	Narrow	Z		904464		904465	904466	904467	904468	
4.3	RePlant™	RP (tri-channel)	T		904374		904375	904376	904377	904378	
3.7	ScrewPlant	3.5 head	B		904340		904341	904342	904343	904344	
4.7	ScrewPlant	4.5 head	C		904345		904346	904347	904348	904349	
IMTEC Corporation®											
3.3	Universal Flare Cylinder	Regular	A				904335	904336	904337	904338	904339
3.75	Universal Self-tapping	Regular	A				904335	904336	904337	904338	904339
3.75	Universal Self-tapping Coated	Regular	A				904335	904336	904337	904338	904339
4.0	Spike Cylinder	Regular	A				904335	904336	904337	904338	904339
4.0	Universal Cylinder	Regular	A				904335	904336	904337	904338	904339
INNOVA											
4.1	ENDOPORE® Ext.	Regular	A				904335	904336	904337	904338	904339

IMPLANT DIAMETER	IMPLANT NAME	PROSTHETIC CONNECTION	GROUP	CUFF HEIGHT						
				0.4mm	1.0mm	1.25mm	2.0mm	3.0mm	4.0mm	5.0mm
4.0	ENTEGRA™ Ext.	Regular	A			904335	904336	904337	904338	904339
Interpore IMZ™										
3.3	Hex Cylinder	Regular	A			904335	904336	904337	904338	904339
3.75	Self-tapping Threaded	Regular	A			904335	904336	904337	904338	904339
4.0	Hex Cylinder	Regular	A			904335	904336	904337	904338	904339
4.0	Self-tapping Threaded	Regular	A			904335	904336	904337	904338	904339
4.25	Hex Cylinder	Regular	A			904335	904336	904337	904338	904339
Keystone (Lifecore)										
3.75	Restore® Self-tapping Screw	Regular	A			904335	904336	904337	904338	904339
4.0	Restore® Self-tapping Screw	Regular	A			904335	904336	904337	904338	904339
3.75	Restore® External Hex Screw	Regular	A			904335	904336	904337	904338	904339
4.0	Restore® External Hex Screw	Regular	A			904335	904336	904337	904338	904339
4.0	Restore® External Hex Cylinder	Regular	A			904335	904336	904337	904338	904339
4.2	Sustain® External Hex Cylinder	Regular	A			904335	904336	904337	904338	904339
3.75	Sustain® External Hex Screw	Regular	A			904335	904336	904337	904338	904339
4.0	Sustain® External Hex Screw	Regular	A			904335	904336	904337	904338	904339
4.2	Sustain® Ext. Hex MC Cylinder	Regular	A			904335	904336	904337	904338	904339
3.3	Stage-1™	4.8 head	S	904264	904265		904266	904267	904268	
4.0	Stage-1™	4.8 head	S	904264	904265		904266	904267	904268	
Minimatic/Stryker										
3.3	External Hex Cylinder	Regular	A			904335	904336	904337	904338	904339
3.75	External Hex Screw	Regular	A			904335	904336	904337	904338	904339
4.0	External Hex Cylinder	Regular	A			904335	904336	904337	904338	904339
4.0	External Hex Screw	Regular	A			904335	904336	904337	904338	904339
4.75	External Hex Screw	Regular	A			904335	904336	904337	904338	904339
5.0	External Hex Cylinder	Regular	A			904335	904336	904337	904338	904339
Nobel Biocare Branemark System®										
3.3	Fixture	Regular	A			904335	904336	904337	904338	904339
3.75	Fixture	Regular	A			904335	904336	904337	904338	904339
4.0	Fixture	Regular	A			904335	904336	904337	904338	904339
5.0	Fixture (Old Version)	Regular	A			904335	904336	904337	904338	904339
3.75	MkII, Self-tapping Fixture	Regular	A			904335	904336	904337	904338	904339
4.0	MkII, Self-tapping Fixture	Regular	A			904335	904336	904337	904338	904339
Nobel Biocare (Steri-Oss®)										
3.8	HL Cylinder	Regular	A			904335	904336	904337	904338	904339
3.8	HL Threaded	Regular	A			904335	904336	904337	904338	904339
4.5	HL Threaded	Regular	A			904335	904336	904337	904338	904339
3.5	Replace Select™	RP (Tri-channel)	Z		904464		904465	904466	904467	904468
4.3	Replace Select™	RP (Tri-channel)	T		904374		904375	904376	904377	904378
5.0	Replace Select™	WP (Tri-channel)	AN		904670		904671	904672	904673	904674
3.5	NobelReplace®	NP (Tri-channel)	Z		904464		904465	904466	904467	904468
4.3	NobelReplace®	RP (Tri-channel)	T		904374		904375	904376	904377	904378
4.0	NobelReplace® Straight	RP (Tri-channel)	T		904374		904375	904376	904377	904378
5.0	NobelReplace®	WP (Tri-channel)	AN		904670		904671	904672	904673	904674
3.5	Nobel Active™	Conical	AY		904680		904681	904682	904683	904684
4.3	Nobel Active™	Conical	AP		904675		904676	904677	904678	904679
5.0	Nobel Active™	Conical	AP		904675		904676	904677	904678	904679
3.5	NobelReplace®	Conical	AY		904680		904681	904682	904683	904684
4.3	NobelReplace®	Conical	AP		904675		904676	904677	904678	904679
5.0	NobelReplace®	Conical	AP		904675		904676	904677	904678	904679
OIC										
3.0	Osteo Standard ST	Regular	A			904335	904336	904337	904338	904339
3.25	Osteo Standard ST	Regular	A			904335	904336	904337	904338	904339

IMPLANT DIAMETER	IMPLANT NAME	PROSTHETIC CONNECTION	GROUP	CUFF HEIGHT						
				0.4mm	1.0mm	1.25mm	2.0mm	3.0mm	4.0mm	5.0mm
3.75	Osteo Standard ST	Regular	A			904335	904336	904337	904338	904339
Sterngold-Implamed®										
3.3	Hex Cylinder	Regular	A			904335	904336	904337	904338	904339
4.0	Hex Cylinder	Regular	A			904335	904336	904337	904338	904339
3.75	Standard Hex Screw	Regular	A			904335	904336	904337	904338	904339
3.75	Self-tapping Hex Screw	Regular	A			904335	904336	904337	904338	904339
3.75	Self-tapping "SST" Hex Screw	Regular	A			904335	904336	904337	904338	904339
4.0	Standard Hex Screw	Regular	A			904335	904336	904337	904338	904339
4.0	Self-tapping Hex Screw	Regular	A			904335	904336	904337	904338	904339
4.0	Self-tapping "SST" Hex Screw	Regular	A			904335	904336	904337	904338	904339
5.0	RP "SST" Hex Screw	Regular	A			904335	904336	904337	904338	904339
3.75	RP Acid Etched	Regular	A			904335	904336	904337	904338	904339
4.0	RP Acid Etched	Regular	A			904335	904336	904337	904338	904339
5.0	RP Acid Etched	Regular	A			904335	904336	904337	904338	904339
3.3	Stern IC	4.8 head	S	904264	904265		904266	904267	904268	
4.1	Stern IC	4.8 head	S	904264	904265		904266	904267	904268	
Straumann										
3.3	ITI TE™	4.8 head	S	904264	904265		904266	904267	904268	
3.3	ITI Std. & Std. Plus	4.8 head	S	904264	904265		904266	904267	904268	
4.1	ITI TE™	4.8 head	S	904264	904265		904266	904267	904268	
4.1	ITI Std. & Std. Plus	4.8 head	S	904264	904265		904266	904267	904268	
4.8	ITI Std. & Std. Plus	4.8 head	S	904264	904265		904266	904267	904268	
4.1	Bone Level	Regular	BD		904469		904470	904471	904472	904474
4.8	Bone Level	Regular	BD		904469		904470	904471	904472	904474
3.3	Bone Level	Narrow	BE		904550		904551	904552	904553	904554
Zimmer (Paragon, Calcitek®, Centerpulse)										
3.5	Bio-Vent® X™	Regular	A			904335	904336	904337	904338	904339
3.75	Swede-Vent™ Conical Neck CST	Regular	A			904335	904336	904337	904338	904339
3.75	Swede-Vent™ Standard	Regular	A			904335	904336	904337	904338	904339
4.0	Swede-Vent™ Standard	Regular	A			904335	904336	904337	904338	904339
4.0	Bio-Vent® X™	Regular	A			904335	904336	904337	904338	904339
3.25	Micro-Vent®	3.5 head	B		904340		904341	904342	904343	904344
3.3	Screw-Vent®	3.5 head	B		904340		904341	904342	904343	904344
3.5	Bio-Vent®	3.5 head	B		904340		904341	904342	904343	904344
3.7	Screw-Vent®	3.5 head	B		904340		904341	904342	904343	904344
3.75	Screw-Vent®	3.5 head	B		904340		904341	904342	904343	904344
4.3	Core-Vent®	3.5 head	B		904340		904341	904342	904343	904344
4.25	Micro-Vent®	4.5 head	C		904345		904346	904347	904348	904349
4.5	Bio-Vent®	4.5 head	C		904345		904346	904347	904348	904349
4.7	Screw-Vent®	4.5 head	C		904345		904346	904347	904348	904349
5.3	Core-Vent®	4.5 head	C		904345		904346	904347	904348	904349
3.7	Tapered Swiss Plus™	4.8 head	S	904264	904265		904266	904267	904268	
4.8	Tapered Swiss Plus™	4.8 head	S	904264	904265		904266	904267	904268	
4.1	Straight Swiss Plus™	4.8 head	S	904264	904265		904266	904267	904268	
4.8	Straight Swiss Plus™	4.8 head	S	904264	904265		904266	904267	904268	
3.75	ThreadLoc™	Regular	A			904335	904336	904337	904338	904339
6.0	Screw-Vent®	5.7 head	BF		904695		904696	904697	904698	904699
Hahn™										
3.5	Tapered Implant		AY							
4.5	Tapered Implant		AY							



Sterngold Dental, LLC
23 Frank Mossberg Drive / Attleboro, MA 02703 USA
1-800-243-9942 or 508-226-5660
www.sterngold.com