

Digital Dentures Built Around **YOU!**

You decide the workflow.

You control your customer communications.

You choose the Digital Denture service options.

You increase and improve overall lab productivity.

You have the **LAB**,
We'll provide the **SUPPORT**.

- Implant Supported Overdenture Integration.
- State-Of-The-Art **Baltic Denture System**[®].
- Merz Dental Denture Teeth Are Pre-Set In **BDLoad**[®].
- Expert Digital Customer Support For Laboratories.










SEE PAGE 2 FOR
DIGITAL WORKFLOW STEPS!

VISIT BDS WEBPAGE
(www.sterngold.com/denture-system)

Simplifying Digital Dentures



Choose a TIER that works for YOU!

	TIER 1 NO DIGITAL EQUIPMENT	TIER 2	TIER 3	TIER 4
STERNGOLD	<ul style="list-style-type: none"> ✓ Scan ✓ Design ✓ Mill ✓ Finish ✓ BDLoad® & BDKEY® Set 	<ul style="list-style-type: none"> ✓ Design ✓ Mill ✓ Finish ✓ BDLoad® & BDKEY® Set 	<ul style="list-style-type: none"> ✓ Mill ✓ Finish ✓ BDLoad® & BDKEY® Set 	<ul style="list-style-type: none"> ✓ BDLoad® & BDKEY® Set
YOU	<p>Select from:</p> <ol style="list-style-type: none"> 1. Completed BDKEY® Set 2. U/L cast with bite registration 3. Impression with the old denture. 	 SCAN Send .STL files	 SCAN  DESIGN*	 SCAN  DESIGN*  MILL  FINISH


*For **DESIGN**, users must send **BDCreator® PLUS** .STL files.

SHORTENED BDS PROCESS

2 appointments needed!



BDKEY® Set
Function and aesthetics in the first appointment.



BDCreator® PLUS
The process-integrated CAD Software.



BDLoad®
The milling blank with occlusive dental arches.



FINISHED BDLoad®
Fully milled, finished and accurate denture.